


## Walkin' Support Sling



### Disclaimer

Please note that this, and all animal support equipment should be used at the advice of a veterinarian or animal rehabilitation specialist. Lifting a dog with this harness can cause problems depending on the health of the dog. In addition, use of a lifting harness can cause chafing. After your initial use, and after each subsequent use, please check the dog carefully for bruises, scrapes, and any chafing.



\*Measure the larger of the dog's abdomen or chest

Size		(A) Abdomen		(B) Chest		(C) Belly	
		inches	cm	inches	cm	inches	cm
Small	from	<b>18</b>	45.72	<b>18</b>	45.72	<b>6</b>	15.24
	up to	<b>22</b>	55.88	<b>22</b>	55.88		
Medium	from	<b>22</b>	55.88	<b>22</b>	55.88	<b>8</b>	20.32
	up to	<b>26</b>	66.04	<b>26</b>	66.04		
Medium/ Large	from	<b>24</b>	60.96	<b>24</b>	60.96	<b>10</b>	25.4
	up to	<b>28</b>	71.12	<b>28</b>	71.12		
Large	from	<b>29</b>	73.66	<b>29</b>	73.66	<b>11</b>	27.94
	up to	<b>35</b>	88.9	<b>35</b>	88.9		
XLarge	from	<b>36</b>	91.44	<b>36</b>	91.44	<b>12</b>	30.48
	up to	<b>44</b>	111.76	<b>44</b>	111.76		