


To measure your dog for the correct size, we will need the girth of one back leg all the way around at the widest part (thigh).  
The Large and XL should be used with the large frame carts



Walkin'<sup>TM</sup>  
Lift Rear Harness



Size	(A) Girth of Leg	
	inches (from-to)	cm (from-to)
XXSmall	1	2.54
	7	17.78
XSmall	7	17.78
	9	22.86
Small	10	25.4
	11	27.94
Medium	12	30.48
	13	33.02
Medium/Large	14	35.56
	15	38.1
Large	16	40.64
	18	45.72
XLarge	19	48.26
	22	55.88

The Walkin' Lift Rear Harness is used to support your dog's hind quarters. This harness is compatible with the Walkin' Wheels Wheelchair, and can be used instead of the Leg Rings. Your Pet can do their business while the harness is on. Handles included. The size of the adjustable handles range from 4.5" to 27".

This Walkin' Lift rear support harness for dogs can also be used as a front harness on dogs that do not have a large chest. It can also be used with the Front Wheel Attachment instead of using the leg rings. Your pet can do it's business in the Walkin' Lift Rear Harness.

**Disclaimer**

Please note that this, and all animal support equipment should be used at the advice of a veterinarian or animal rehabilitation specialist. Lifting a dog with this harness can cause problems depending on the health of the dog. In addition, use of a lifting harness can cause chafing. After your initial use, and after each subsequent use, please check the dog carefully for bruises, scrapes, and any chafing.